

Building experiments in



By Johnny Lau

Modified on 26 Mar, 2019

PsychoPy – Exercise

Stroop task

Task1 (~15mins)

- Familiarising with the software
- Basic control



Task 2 (~25mins)

- Exploring other features
- Coding

Task 2

- Exploring other components (Image)
- Simple coding
- Introduce breaks in the experiment


Illustration of the experimental task




Preparation

To complete this task, you can either modify your programme from Task1 or start fresh by using the folder '**Ex2_stroop_task-imgver**'

In the folder '**Ex2_stroop_task-imgver**', there is a subfolder 'img' which contains all the image files you need to complete Task 2. There are 8 images altogether.



Conditions file



Similar to task 1, you need a conditions file so that the programme would know what stimulus items you are actually presenting

- Now, open the '**conditions2**' excel file
- You probably notice the first 4 columns contain exactly the same information as in the conditions file in Task 1 (i.e. **word, colour, congruent, corrAns**)
- In addition, in this '**conditions2**' file, we need an extra column '**img_path**' to record where the image files are located

	A	B	C	D	E	F	G	H
1	word	colour	congruent	corrAns	img_path			
2	red	red	1	left	img/red-red.png			
3	red	blue	0	right	img/red-blue.png			
4	blue	red	0	right				
5	blue	blue	1	left				
6	white	black	0	right				
7	black	white	0	right				
8	white	white	1	left				
9	black	black	1	left				
10								
11								
12								
13								
14								
15								

You need to complete this column with the corresponding image paths accordingly

A bit more info. about specifying the path of an image

	A	B	C	D	E	F	G	H
1	word	colour	congruent	corrAns	img_path			
2	red	red	1	left	img/red-red.png			
3	red	blue	0	right	img/red-blue.png			
...	..	.	-	...				

For example, for the image **red-blue.jpg** (which is the word 'RED' coloured in blue), its path should be written like this in the conditions file:

img/red-blue.png


('img' folder is on the same level as your programme file and you will need to specify where your image is located ONLY from this level)

*You don't need to write out its full path (e.g. C:\Username\Desktop\WS3_sempriming\...) because PsychoPy will automatically fill this out for you


Now, open the **'try2_stroop.psyexp'** programme file (builder view) in the folder **'Ex2_stroop_task-imgver'**.

The main trials

Creating the main trials:



How this would be represented in PsychoPy:




The main trials: Target image

Update the fields in Image Properties as follows

The screenshot shows the 'Image Properties' dialog box in a software interface. The dialog box has the following fields and settings:

- Name:** image
- Start:** time (s) [blank] (An annotation points to this field: "Leave this blank for now")
- Expected start (s):** [blank]
- Stop:** duration (s) 3
- Expected duration (s):** [blank]
- Image:** \$img_path (An annotation points to this field: "Image: \$path")
- Position [x,y]:** (0,0)
- Size [w,h]:** [6,2.5] (An annotation points to this field: "Play around these measures to achieve the best dimension for presenting the image")
- Orientation:** 0
- Opacity:** 1
- Units:** cm
- Repeat:** set every repeat (An annotation points to this dropdown: "**choose among the drop-down options: 'set every repeat' (this changes trial after trial)")
- Constant:** constant

The background shows a software interface with a menu bar (Edit, Tools, View, Experiment) and a toolbar. The main area displays a trial sequence: Instructions, trial, breaktime_2, and EndingMessage (6.00s). The trial sequence is labeled 'trials_loop'.



This time we want the fixation cross to be presented with a variable length of duration each time (i.e. different duration in a different trial)

With this line, you create a variable `fix_dur` which varies in value between 1 and 2 each time

This can be achieved by writing some simple codes in a Code component

The screenshot shows the PsychoPy software interface with the following elements:

- Code Component:** A component named `code_2` is circled in red. It contains the following code:


```
y=sin(x)
x=random()
f='tst.jpg'
```
- Code Editor:** A window titled "Begin Routine *" is open, showing the following code:


```
#generate a new variable each time this routine is loaded
#google what the randint function in Python does
fix_dur=randint(1,3)

#this line saves the variable generated in the output file
thisExp.addData('fix_duration', fix_dur)
```

 A red arrow points from the text "With this line, you create a variable `fix_dur` which varies in value between 1 and 2 each time" to the line `fix_dur=randint(1,3)`.
- Routine Loop:** At the bottom, a flowchart shows the sequence of components: `Instructions` → `trial` → `breaktime` → `EndingMessage (6.00s)`. A `trials_loop` box is positioned below the `trial` and `breaktime` components, with arrows indicating a loop between them.
- Stimuli and Responses:** On the right side, there are panels for "Stimuli" and "Responses". The "Responses" panel has a "Custom" component circled in red, containing the code:


```
y=sin(x)
x=random()
f='tst.jpg'
```

fix_cross Properties

Basic

Name: fix_cross

Start: time (s) 0.0

Stop: duration (s) fix_dur

Color: black

Font: Arial

Letter height: 0.12

Position [x,y]: [0, 0]

Text

Help

Now, define the duration as variable 'fix_dur'

image Properties

Basic

Name: image

Start: time (s) fix_dur

Stop: duration (s) 3

Image: \$img_path

Position [x,y]: (0,0)

Size [w,h]: [6,2.5]

Orientation: 0

Opacity: 1

Units: cm

Help OK Cancel

trial_resp Properties

Basic

Name: trial_resp

Start: time (s) fix_dur

Stop: duration (s) 3

Force end of Routine: ☒

Allowed keys: 'left', 'right'

Store: last key

Store correct: ☒

Correct answer: \$corrANS

Discard previous: ☒

sync RT with screen: ☒

Help OK Cancel

And the response key starts only after the fix_cross finishes i.e. after the time defined by variable 'fix_dur'

The main trials: Loops

Just similar to the previous example, Task 1

Between trials: give me a break!

Inserting breaks between trials




- Introducing breaks among trials = splitting the trials into various parts
- You will need to create a new routine and perhaps name it '**breaktime**'. This 'breaktime' routine should also be wrapped around inside the trials loop.


- Click Insert Routine > (new)
- Type '**breaktime**'
- Place it after **trial** routine and make sure it is wrapped within the loop

Inside the breaktime routine, you will need these components

Breaks




1. First, insert a **Text** component
Update the Properties as shown on the left





2. Then, use a Keyboard component so that the participant can respond to restart the experiment at their own pace
Update the Properties as shown on the right

Breaks




3. Finally, you need to use a **'code'** component
- On **'Each Frame'** tab, type these codes:

*Make sure this part here is exactly the same as the name of the loop (i.e. **trials_loop**) you've created*



We will have in total 80 trials (i.e. 8 items x 10 repetitions).


Here, we want to introduce seven breaks (after each 10 trials)




At the end, your **'breaktime'** routine panel should look like this.

Also, it's important to place the **code** component at the very bottom among all

Now test your programme!



Other resources



2. Official website:

- <http://www.psychopy.org/>

3. Forum:

- <https://discourse.psychopy.org/>

4. Google!

PsychoPy Online

<https://www.psychopy.org/online/fromBuilder.html>

- **Sync your files with Pavlovia.org**

<https://www.pavlovia.org/>